

(3mm aperture, 3000nm to 2000nm, DC-MHz)



**DATASHEET** 

Return to the Webpage \( \int\_0 \)




### **Features**

- High Performance
- Compact Package
- Easy integration
- Customize Available
- Low cost

## **Applications**

- Laser Modulation
- Holography
- Metal cutting/engraving
- Microfabrication

The FEOM Series Free-space Electro-Optic Modulator, supplied with a matched driver, modifies the polarization or amplitude of free-space laser beams across a wide wavelength range. It uses a compensated electro-optic crystal pair with specialized material compositions for high-power handling and short-wavelength stability. To lower the driving voltage, the modulator is offered in two configurations: single-pair and dual-pair crystal lengths, as shown in the accompanying graph. The device operates with alternating voltage to prevent internal charge buildup, while the matched driver supports modulation from DC to MHz, with depth limited by the driver's output power. Optional accessories include polarization cubes for high-extinction-ratio amplitude modulation and isolators to block back reflections. A polarized input beam is required for amplitude modulation.



### **Specifications**

Parameter	Min	Typical	Max	Unit		
Wavelength Range	W1	400		600		
	W2	600		900		
	W3	900		1250	nm	
	W4	1250		1650		
Clear Aperture	3			mm		
Halfwave Voltage, non-resonant	80		600	٧		
Extinction Ratio [1]	10	20	30	dB		
Input impedance, resonant		50		ohms		
Input capacitance, non-resonant		14		pF		
Max Optical Power Density	532nm		2	10 [2]	W	
	1064nm		5	20 [2]	W	
Temperature		-20		50	°C	

#### lotes:

- [1]: Related to the polarizer and beam quality
- [2]: High power version available

**Note:** The specifications provided are for general applications with a cost-effective approach. If you need to narrow or expand the tolerance, coverage, limit, or qualifications, please [click this <u>link</u>]:

Legal notices: All product information is believed to be accurate and is subject to change without notice. Information contained herein shall legally bind Agiltron only if it is specifically incorporated into the terms and conditions of a sales agreement. Some specific combinations of options may not be available. The user assumes all risks and liability whatsoever in connection with the use of a product or its application.

Rev 07/11/25

© Photonwares Corporation

P +1 781-935-1200

sales@photonwares.com

www.agiltron.com




(3mm aperture, 3000nm to 2000nm, DC-MHz)



### **DATASHEET**

## **Mechanical Drawing (mm) 1X**



<sup>\*</sup>Product dimensions may change without notice. This is sometimes required for non-standard specifications.

## **Operation Instructions**

#### 1. Input Alignment

Direct the laser beam through the device, ensuring proper alignment along the optical axis without hitting the crystal walls.

### 2. Input Polarizer Setup

Place a vertical polarizer at the input side. Adjust the input laser polarization direction to achieve maximum output intensity.

#### 3. Output Polarizer Setup

Place a horizontal polarizer at the output side of the device. Carefully adjust the orientation of the output polarizer to achieve minimum transmitted output intensity (extinction condition).

### 4. Apply Control Voltage

Gradually apply voltage to the device and observe changes in the output intensity. At the operating voltage (Vp), the output intensity will reach its maximum.




(3mm aperture, 3000nm to 2000nm, DC-MHz)




DATASHEET

## **Mechanical Drawing (mm) 2X**



<sup>\*</sup>Product dimensions may change without notice. This is sometimes required for non-standard specifications.







(3mm aperture, 3000nm to 2000nm, DC-MHz)




**DATASHEET** 

## **Modulator Half-Wave Voltage**





## **Amplitude Electro-Optic Crystal Configuration (yellow indicates electrode)**



## **Typical Resonance Response (sine wave)**





(3mm aperture, 3000nm to 2000nm, DC-MHz)



## **DATASHEET**

## **Ordering Information**

Prefix	Туре	Wavelength	Optical Power	Crystal Length	Input Cube **	Output Cube **	Driver	Isolator
<b>EOMN</b> - FEOM*	Amplitude = 5 Phase = 6	400~600 nm = 05 600~900 nm = 07 900~1250 nm = 09 1250~1650 nm = 14 320~500 nm = 03 1920-2400 nm = 20	Regular = 1 High Power = 2	1x = 1 2x = 2 3x = 3	No = 1 Polacore = 3 PBS = 4 Glan-Thompson = 5	No = 1 Polacore = 3 PBS = 4 Glan-Thompson =5	Non = 1 Yes = 2	None = 1 One Stage = 2 Two Stage = 3

<sup>\*</sup> Old part number

## **Polarizer's prices:**

Polacore	\$256		
PBS	\$240		
Glan Thompson	\$485		

### **Polarizer**







<sup>\*\*</sup> Polacore - CW 10W/cm<sup>2</sup> PBS - CW 15W/cm<sup>2</sup>  $\textbf{Glan-Thompson} - \text{CW } 2\text{kW/cm}^2$